
MAWe

by GADNIC

SACALECHE
ELÉCTRICO
SILKE1

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BEBE0001

Índice

Español.....	3
Português.....	10
English.....	17

Características funcionales

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Functionalities

Product Features

1. Automatic simulation: Controlled by a microchip program, it mimics babies' intermittent breastfeeding.
2. In galactagogue mode, the intensity can be adjusted in 9 levels. In massage mode, the intensity can be adjusted in 9 levels to effectively massage the breast. In breast suction mode, the intensity can be adjusted in 9 levels, mimicking the speed of the baby's suction rhythm and the size of the baby's suction volume.


• Product Features

1. Cylinder design with gas-liquid separation structure, no milk backflow when sucking.
2. Removable liquid silicone horn mimics the principle of ergonomics. It can suck and massage, tighten the breast, and stimulate breastfeeding while sucking.
3. Touch-sensitive control panel design, display interface, clear operation process.

Introduction of the pieces

Piece name	Quantity	Material	Piece name	Quantity	Material
Dirt cover	2	PP	Silicone speaker	2	Silica gel
Horn mouth	2	PP	Silicone cylinder	2	Silica gel
Three-way	2	PP	Cylinder head	2	PP
Assembly	4	PP/TPR	Flexible tube	2	Silica gel
Cable	1	/	Duckbill valve	2	Silica gel
Main motor	1	/	Pacifier kit	2	Silica gel+PP
Wide-mouth bottle of 180ml	2	PP			

Next, the parts for the single-bottle milk extractor with only 1 piece

Malfunction analysis

The following problems may occur during milk extraction. Please refer to the provided method.

Failure	Analysis of causes	Processing method
Weak suction	The connection has air leaks (duckbill valve, silicone nozzle, three-way main motor, and straight head). The air outlet route of the main motor is blocked, such as hose folding or debris blockage in the hose airway.	Check the connection parts, reinstall and clean them.
No suction	Duckbill valve not installed or deformed. Improper use of plastic components causing deformations and damage.	Install new duckbill valves or repair damaged accessories.
No suction	The adapter is not connected to the socket or the main unit socket.	Reconnect the adapter before use.
No power	Check for milk or water entering the cylinder or hose to verify if the cylinder is damaged or not properly installed.	Stop immediately, clean or replace the cylinder and reinstall it.

Maintenance

- Do not disassemble or dismantle the product.
- The host power supply requires a 5V DC power adapter or a specified mobile phone power adapter. Do not use other non-standard power adapters to power it.
- Do not place near high temperature sources or fire during use, and avoid placing the main unit on a heater or fire source to prevent combustion and explosion.
- If there is a failure in the mainframe, do not disassemble or attempt to repair it.
- It must be sent to our company's designated maintenance point for repair or replacement by professional personnel.
- Before use, accessories directly in contact with breast milk (bottle, silicone horn, three-way, cylinder, duckbill valve, bottle cap) must be cleaned and disinfected in advance. They can be sterilized with steam or hot water, and the disinfection time should not exceed 1 minute. The host should not be submerged in water.

Screen display

- When the working mode is prolactin,  the light is on; when the working mode is "massage", the light  is on; when it is turned on , it is in working mode; when the "x" is on, it means off.
- During the "prolactin" mode, the intensity can be divided into 9 levels, and the screen will also display the corresponding number. In "Massage" mode, there are 9 levels, and the screen will display the corresponding number.
- In suction mode, there are nine levels, and the screen displays the corresponding numbers.

Regulation of suction strength frequency

- In "Feeding" mode, pressing + and - can adjust the milk strength, a total of 9 options. In
- "Massage" mode, pressing + and - can adjust the prolactin intensity, a total of 9 levels. In "Suction" mode, by pressing + and - you can adjust the suction strength, a total of 9 levels. During operation, pressing the M button can freely change the working mode.
- When switching to "Prolactin" mode, the default value is the last used level; when switching to "Suction" mode, the default value is the last used level.

Suspension function

- Press the pause mode, after this, feel free to press the button to resume.

Memory function


- Start in prolactin mode, massage the suction mode, adjust the required gear for more than 3 seconds, it will automatically remember. The next startup will resume the previous setting.

Power


- In the operating state, the power adapter is connected and powered by the power source. When the power adapter is not connected, the host will shut down.

Tip: Before each use, it must be sterilized beforehand with breast milk that comes into direct contact with the components: bottle, silicone horn, three-way silicone cylinder, duckbill valve.


Bottle cap, and so on


1. Place the silicone horn and bell mouth in sequence, then place the dust cover.

2. Insert the set of silicone horn and bell mouth into the three-ways to ensure the sleeve is tight.
(Note: when removing, one hand holds the three-ways, the other hand holds the bell mouth assembly upwards and gently pulls out).


3. Place the cylinder and cylinder head into three-ways and tighten the cylinder head.

4. Place the duckbill valve on the lower outlet of the three-ways valve to ensure no leaks.

5. Screw the three-ways assembly onto the bottle.


6. Insert the direct hose assembly into the cylinder head and the relative position of the main motor hole to ensure no air leaks and it can be used.

Note: The pass-through assembly must be opened before insertion into the main airway hole; Do not start the main motor work without opening the dust cover, otherwise, it will affect the main motor's service life.

Ready to use

- Before using, check and make sure that the sterilized parts (parts in direct contact with breast milk) are clean and free of water, and that the absorber is firmly connected to the main structure.
- Clean, compress, and massage the breasts and areola with warm, damp towels before suctioning. Relax your body to alleviate the discomfort of engorged breasts.

Panel operation instructions

Startup mode switch

- Touch "⏻" for 2 seconds and enter the default prolactin mode. After startup, you can
- proceed with prolactin, massage, and suction mode.

Rev.1