# GADNIC

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
## **Blood Pressure Monitor**

**Wrist Type**

Thank you for purchasing the Blood Pressure Monitor. The unit has been constructed using reliable circuitry and durable materials. Used properly, this unit will provide years of satisfactory use.

This device is intended for non-invasive measuring of an adult individual's systolic, diastolic blood pressure and heart rate using the oscillometric method. The device is not intended for use on infants and children. The device is designed for home or clinical use. All functions can be used safely and values can be read out in one LCD DISPLAY. Measurement position is on adult wrist only.

Please read this manual thoroughly before using the unit. Please retain this manual for future reference. For specific information about your blood pressure, please CONSULT YOUR DOCTOR. The PATIENT is an intended OPERATOR.

To avoid risk and damage follow all warning precautions.  
Operate unit only as intended. Read all instructions prior to use.

## WARNING SIGNS AND SYMBOLS USED

	Caution
	Mandatory
	Prohibited
	Type BF Equipment
	Instructions For Use MUST be Consulted
	Serial Number
	Discard the used product to the recycling collection point according to local regulations
	The product conforms to the requirements of the Regulation (EU) 2017/745 MDR on medical devices
	Manufacturer
	Authorised Representative in the European Community
	Keep off Sunlight
	Manufacturing Date
	Medical Device

## Caution

Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.

Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting your physician first.

Product is designed for its intended use only.  
Do not misuse in any way.

Product is not intended for infants or individuals who cannot express their intentions.

Do not disassemble or attempt to repair.

Do not use cell phones and other devices, which generate strong electrical or electromagnetic fields, near the device, as they may cause incorrect readings and interference or become interference source to the device.

## Battery Precautions

Do not mix new and old batteries simultaneously.

Replace batteries when Low Battery Indicator "  " appears on screen.

Be sure battery polarity is correct.

Do not mix battery types. Long-life alkaline batteries are recommended.

Remove batteries from device when not in operation for more than 3 months.

Dispose batteries properly; observe local laws and regulations.


## Important Instructions Before Use

1. Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
2. Contact your physician if test results regularly indicate abnormal readings.
3. If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.
4. Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.
5. For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
6. People suffering from vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant should consult their physician before measuring their blood pressure themselves. Different values may be obtained due to their condition.
7. People suffering from arrhythmias such as atrial or ventricular premature beats or atrial fibrillation only use this blood pressure monitor in consultation with your doctor. In certain cases oscillometric measurement method can produce incorrect readings.


- 8.Too frequent measurements can cause injury to the patient due to blood flow interference.
- 9.The cuff should not be applied over a wound as this can cause further injury.
- 10.**DO NOT** attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.
- 11.The cuff should not be placed on the arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.
- 12.Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
- 13.A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.
- 14.Check that operation of the unit does not result in prolonged impairment of the circulation of the patient.
15. Product is designed for its intended use only. Do not misuse in anyway.
16. Product is designed for its intended use only. Do not misuse in any way.
17. Product is not intended for infants or individuals who cannot express their intentions.
18. Prolonged over-inflation of the bladder may cause ecchymoma of your arm.

18. Do not disassemble the unit or wrist cuff. Do not attempt to repair.
19. Use only the approved wrist cuff for this unit. Use of other wristcuffs may result in incorrect measurement results.
20. The system might produce incorrect readings if stored or used out-side the manufacturer's specified temperature and humidity ranges.  
Make sure to store the blood pressure monitor, children, pets and pests are outside of accessible range.
21. Do not use the device near strong electrical or electromagnetic fields generated by cell phones or other devices, they may cause incorrect readings and interference or become interference source to the device.
22. Do not mix new and old batteries simultaneously.
23. Replace batteries when Low Battery Indicator “ ” appears on screen.  
Replace both batteries at the same time.
24. Do not mix battery types. Long-life alkaline batteries are recommended.
25. Remove batteries from device when not in operation for more than 3 months.
26. Do not insert the batteries with their polarities incorrectly aligned.
27. Dispose batteries properly; observe local laws and regulations.
28.  Advising operator that Instruction manual/ Booklet must be consulted.
29. The PC with connection to the device with USB shall meet the requirements of standard IEC 60601-1 or IEC 60950-1.
30. Do not use the device during transport vehicles for influencing measurement accuracy,such as patient transport in an ambulance or helicopter.
31. Contains small parts that may cause a chocking hazard if swallowed by infants.

## Monitor Unit


## Display


## Contents

1. Monitor Unit

2. Owner's Manual

3. Plastic Storage Case

# Important Testing Guidelines

126

1. Avoid eating, exercising, and bathing for 30 minutes prior to testing.
2. Sit in a calm environment for at least 5 minutes prior to testing.
3. Do not stand while testing. Sit in a relaxed position while keeping your wrist level with your heart.
4. Avoid speaking or moving body parts while testing.
5. While testing, avoid strong electromagnetic interference such as microwave ovens and cell phones.
6. Wait 3 minutes or longer before re-testing.
7. Try to measure your blood pressure at the same time each day for consistency.
8. Test comparisons should only be made when monitor is used on the same wrist , in the same position, and at the same time of day.
9. This blood pressure monitor is not recommended for people with severe arrhythmia.
10. Do not use this blood pressure monitor if the device is damaged.

# Quick Start

127

1. Install batteries. (See Figure A )


Figure A

2. Remove clothing from the wrist area. (See Figure B)
3. Rest for several minutes prior to testing. Wrap cuff around left wrist . (See Figure C)


Figure B


Figure C

4. Sit in a comfortable position and place wrist level with heart.  
(See Figure D)
5. Press "■" button to start testing. (See Figure E)

Figure D



Figure E

## Battery Installation

Slide battery cover off as indicated by arrow.

Install 2 new AAA alkaline batteries according to polarity.

Close battery cover.

Note:

- 1) Replace batteries when Low Battery Indicator "  " appears on screen.
- 2) Batteries should be removed from device when not in operation for an extended period of time.


## System Settings

With power off, press and hold "  " button to actuate system setting .The

Memory Group icon flashes.

### 1. Select memory Group

While in the System Setting mode you may accumulate test results into 2 different groups. This allows 2 different users to save individual test results. Press "  " button to choose a group setting. The test results will automatically store in each selected group.



### 2. Time/Date setting

Press "  "button again to enter set the year mode. The corresponding parameters can be set by adjusting the"  " button.

Press "  " button again to enter set the month mode. The User can continue setting the day, hour ,minute and 12/24hours in the same way. Every time the "  " button is pressed , it will lock in your selection and The corresponding parameters can be set by adjusting the "  "button. ( month,day,hour, minute,12/24 hours)


# Unit Operation

131

## 3. Time Format Setting.

Press "  " button again to set the time format mode. Set the time format by adjusting the "  " button. EU means European Time US means U.S Time.

## 4. Voice Setting

Press "  " button to enter voice setting mode. Set the voice format ON or OFF by pressing the "  " button.

## 5. Volume Settings

Press "  " button to enter volume setting mode.


Set the voice volume by adjusting the "  " button.

Smaller "  " is for lower volume. There are six volume levels.

## 6. Settings of cuff keeping level with heart

Press "  " button to enter Settings of cuff keeping level with heart .Set the Function of cuff keeping level with heart ON or OFF by pressing the "  " button.

## 7. Save Settings


While in any setting mode, press "  " button to turn the unit off. All information will be saved.

Note: Unit will automatically save all information and shut off if left idle for 3 minutes.

## Applying The Wrist Monitor

Do not apply over clothing. If wearing a long sleeved shirt, be sure to roll sleeve back to forearm.

Apply monitor to wrist as illustrated. Tighten cuff firmly as not to wiggle.

Do not stand while testing. Sit in a comfortable position with back supported, feet flat on the floor with legs uncrossed. Place middle of the cuff at the level of the right atrium of the heart.

## Testing

### 1. Power On


Press “” button to turn on the unit . The LCD screen will appear full display for one second as unit performs a quick diagnosis. The unit is ready for testing.


Note: Unit will not function if residual air from previous testing is present in cuff. The LCD will flash "↓" until pressure is stabilized.

## 2. Make sure the wristband is level with the heart

When the wristband does not keep level with the heart, the LCD screen displays "  " and flashes, as shown in the following interface (the buzzer "di di" prompts the user current location error )


When the wrist strap and the heart remain level, enter the zero grasp, the LCD screen will display as follows interface:

### 3. Testing

After cuff inflation, air will slowly rise as indicated by the corresponding cuff pressure value. A flashing "  " will appear simultaneously on screen signaling heart beat detection.

Note: Remain relaxed during testing. Avoid speaking or moving body parts.

## 4. Result Display

The screen will display measurements for systolic and diastolic blood pressure. An indicator representing the current measurement will appear next to the corresponding WHO Classification.

## Irregular Heartbeat Indicator

If the monitor detects an irregular heart rhythm two or more times during the measuring process, the Irregular Heartbeat Symbol " (Heart) " appears on screen along with measurement results. Irregular heartbeat rhythm is defined as rhythm that is either 25% slower or faster than the average rhythm detected while measuring systolic blood pressure and diastolic blood pressure. Consult your physician if the Irregular Heartbeat Symbol " (Heart) " frequently appears with your test results.

## Power Off

The "  " button can be pressed to turn off the unit in any mode.


The unit can turn off the power itself about 3 minutes no operation in any mode.

Safety Precaution: If pressure in cuff becomes too extreme while testing, press the "  " button to turn power off.


The cuff pressure will rapidly dissipate once the unit is off.

## Memory Check and Last 3 Tests Average

With power off, press the "■" button to activate screen display. After the unit performs a self diagnosis, the screen will display the average test results from the last 3 readings of the last group used . The " **AVG**" symbol will appear along with the corresponding WHO Blood Pressure Indicator. To check the average results from other groups, select the desired group first prior to activating the "■" button in the off position (See "Select Memory Group" on page 16).

Press the "➡" button again, you may check past test results. Upon activating test results, you can press the "➡" button to scroll through all test results stored in memory. The LCD will display the last memory as NO: 01 reading. ("➡" button means turn forward,)


Note: Past test results will only be displayed from the most recently used memory group. To check past test results in other memory groups, you must Press corresponding button or select the desired group and then turn monitor off .

## Memory Deletion

Memory for a selected group may be deleted while in Memory

Check mode. Press and hold the "  "button 3 second for approximately 3 seconds to delete all memory records from the selected group with voice broadcast " Memory Clear ". And then transfer into testing mode. Press the "  " button to turn the unit off.

Note: Memory cannot be recovered once it has been deleted.

## Low Battery Indicator

The unit will broadcast " Low Battery"when battery life is depleting and unable to inflate cuff for testing. The "  " appears simultaneously for approximately 5 seconds prior to shutting off.

Replace batteries at this time. No memory loss will occur throughout this process.


## Arm Shake Indicator

If there is arm movement during the measurement, the "  " icon may flash. Indicates that the measurement results may be inaccurate, and the situation will be recorded at the end of the measurement as a reminder.

## Static Pressure Measurement

In the power down state, press and hold the “” button 3 second, and then install the batteries. until the LCD screen is full, release the “” button. When the LCD screen displays the double zero, the blood pressure meter is in static state. Software version is displayed 10 is a software version in the figure.

Note: Only Service personnel permitted to access to this mode,  
the mode unavailable in normal use.

# Unit Operation

## Troubleshooting

Abnormal phenomenon	Cause analysis	Processing method
Abnormal sphygmomanometer	The armband is tied too tight or too loose; Or the arm strap is tied incorrectly;	Roll the armband correctly
	Move the arm during measurement or Electronic sphygmomanometer	Stay quiet, keep your arm steady, and do not move the monitor
	Speaking, nervous or emotional during measurement	Instead of talking, take deep breaths to calm your mood and relax your body
	Incorrect measurement posture	Adjust posture, see "Blood pressure gauge Wearing"
	There is interference in charging process or improper operation in measuring process	See operation Instructions.

The following table shows the error signs that may occur during measurement, possible causes and handling methods. Please measure again using the correct method

Error display	The cause of the problem	The solution
Er1	Can't detect high and low pressure	Please fasten the cuff before measuring
Er2	Cuff too loose or loose	Please fasten the cuff before measuring
Er3	Improper compression caused by arm or body movement	Hold the arm or body still and measure again
Er4	The pressure exceeds 300mmHg	Please fasten the cuff before measuring
Er5	The pressure exceeds 15mmHg for 3 minutes	Check whether the cuff is knotted or the vent valve is blocked. If the problem persists, contact the manufacturer
Er6	Blood pressure measurements were out of range	Please tighten and measure and measure again If you cannot solve the problem, please contact the manufacturer
	Battery dead	Replace the battery or connect the power adapter (if any).


Note: If you cannot solve the abnormal situation by yourself, you can consult the manufacturer or the manufacturer's designated unit by phone. It is forbidden to disassemble and repair without permission. If necessary, professional maintenance personnel can ask the manufacturer for the list of components and circuit schematic diagram.

## Blood Pressure

Blood pressure is the force of blood pushing against the walls of arteries. It is typically measured in millimeters of mercury (mmHg.) Systolic blood pressure is the maximum force exerted against blood vessel walls each time the heart beats. Diastolic blood pressure is the force exerted on blood vessels when the heart is resting between beats.

An individual's blood pressure frequently changes throughout the course of a day. Excitement and tension can cause blood pressure to rise, while drinking alcohol and bathing can lower blood pressure. Certain hormones like adrenaline (which your body releases under stress) can cause blood vessels to constrict, leading to a rise in blood pressure.



If these measuring numbers become too high, it means the heart is working harder than it should.

## WHO Blood Pressure Classification Indicator

The product is equipped with a classification indicator based on established guidelines from the World Health Organization.


The chart below (color coded on monitor unit) indicates test results.

↳: Blood Pressure Classification Indicator

## Health Reminder

Hypertension is a dangerous disease that can affect the quality of life. It can lead to a lot of problems including heart failure, kidney failure, and cerebral hemorrhaging. By maintaining a healthy lifestyle and visiting your physician on a regular basis, hypertension and relative diseases are much easier to control when diagnosed in the early stages.

Note: Do not be alarmed if an abnormal reading occurs.

A better indication of an individual's blood pressure occurs after 2-3 readings are taken at the same time each day over an extended period of time. Consult your physician if test results remain abnormal.

**Q:** What is the difference between measuring blood pressure at home or at a professional healthcare clinic?

**A:** Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.

Note: Abnormal test results may be caused by:

1. Improper cuff placement

Make sure cuff is snug-not too tight or too loose.

2. Improper body position

Make sure to keep your body in an upright position.

3. Feeling anxious or nervous

Take 2-3 deep breaths, wait a few minutes and resume testing.

**Q:** What causes different readings?

**A:** Blood pressure varies throughout the course of a day. Many factors including diet, stress, cuff placement, etc. may affect an individual's blood pressure.


**Q:** Should I apply the cuff to the left or right wrist? What is the difference?

**A:** Either wrist can be used when testing, however, when comparing results, the same wrist should be used. Testing on your left wrist may provide more accurate results as it is located closer to your heart.

**Q:** What is the best time of day for testing?

**A:** Morning time or any time you feel relaxed and stress free.

1. Avoid dropping, slamming, or throwing the unit.

2. Avoid extreme temperatures. Do not use outdoors.


- When cleaning the unit, use a soft fabric and lightly wipe with mild detergent. Use a damp cloth to remove dirt and excess detergent.


- Cuff Cleaning: Do not soak cuff in water! Apply a small amount of rubbing alcohol to a soft cloth to clean cuff's surface. Use a damp cloth (water-based) to wipe clean. Allow cuff to dry naturally at room temperature. The cuff must be cleaned and disinfected before use between different users.


- Do not use petrol, thinners or similar solvents.

- Remove batteries when not in operation for an extended period of time.

7. Do not disassemble product.

8. It is recommended the performance should be checked every 2 years.
9. Expected service life: Approximately three years at 10 tests per day.
10. No service and maintenance while it is in use and maintenance only be performed by service personnel. Service and maintenance require parts, repair, technical support will be provided.

# Specifications

157

Product Description	Wrist-type Fully Automatic Blood Pressure Monitor	
Display	LCD Digital Display Size:45mm x 30.5mm(1.77" x 1.20")	
Measurement Method	Oscillometric Method	
Measurement Range	Systolic Pressure	60mmHg~260mmHg
	Diastolic Pressure	40mmHg~200mmHg
	Pressure	0mmHg~299mmHg
	Pressure	±3mmHg
	Pulse	30 ~ 180 Beats/Minute
	Pulse	±5%
Pressurization	Automatic Pressurization	
Memory	2X150 Memories in Tow Groups with Date and Time	
Function	Irregular Heartbeat Detection	
	WHO Classification Indicator	
	Last 3 Results Average	
	Low Battery Detection	
	Automatic Power-Off	

# Specifications

158

Function	Voice
	Backlight
Power Source	2 Alkaline Batteries Size AAA
Battery Life	Approximately 2 months at 3 tests per day
Unit Weight	Approx. 91g (3.21 oz) (Excluding Battery)
Unit Dimensions	Approx. 86.5mm×64.7mm×26.6mm(L x W x H) (3.41" x 2.55"x 1.05")
Cuff Circumference	Fits wrist circumference 13.5-21.5 cm(5.3"-8.5")
Operating Environment	Temperature 10°C ~ 40°C (50°F~104°F)
	Humidity 15%~93%RH
	Pressure 800hPa~1060hPa
Storage Environment	Temperature -25 °C~55 °C (-13°F~131°F)
	Humidity ≤93% RH
Transport Environment	Temperature -25 °C~55 °C (-13°F~131°F)
	Humidity ≤93% RH
Ingress Protection Rating	IP 22
Classification	Internal Powered Equipment Type BF 

Specifications are subject to change without notice.

This Blood Pressure Monitor complies with the European regulations and bears the CE mark “CE 0123”. This blood pressure monitor also complies with mainly following standards

(included but not limited):

Safety standard:

EN 60601-1 Medical electrical equipment part 1: General requirements for safety EMC standard:

EN 60601-1-2 Medical Electrical Equipment -- Part 1-2: General Requirements For Basic Safety And Essential

Performance -- Collateral Standard:

Electromagnetic Disturbances –Requirements And Tests.

Performance -- Collateral Standard:

Electromagnetic Disturbances –Requirements And Tests.

Performance standards:

IEC80601-2-30, Medical electrical equipment – Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers.

EN 1060-3 Non-invasive sphygmomanometers - Supplementary requirements for electromechanical blood pressure measuring systems.

ISO 81060-2, non-invasive sphygmomanometers - part 2: clinical validation of automated measurement type.

The Blood Pressure Monitor is guaranteed from the date of purchase. If the Blood Pressure Monitor does not function properly due to defective components or poor workmanship, we will find a possible solution. The warranty does not cover damages to your Blood Pressure Monitor due to improper handling. Please contact local retailer for details.

# Electromagnetic Compatibility Information

162

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

**Table 1**

Guidance and declaration of manufacturer-electromagnetic emissions	
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.	
Emissions test	Compliance
Radiated emission CISPR 11	Group 1, ClassB

**Table 2**

Guidance and declaration of manufacturer-electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
IMMUNITY test	IEC 60601 test level	Compliance level
Radiated RF EM fields IEC 61000-4-3	3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz (See table 3)	3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz (See table 3)
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz
Electrostatic-discharge (ESD) IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Proximity magnetic field IEC 61000-4-39	See Table 5	See Table 5

## **Electromagnetic Compatibility Information**

164

Table 3

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device		
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = [\frac{3.5}{E_1}] \sqrt{P}$	800 MHz to 2.7 GHz $d = [\frac{7}{E_1}] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**Table 5**

Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields

Test frequency	Modulation	IMMUNITY TEST LEVEL(A/m)
30 kHz a)	CW	8
134,2kHz	Pulse Modulation $\approx$ 2,1kHz	65 c)
13,56MHz	Pulse Modulation $\approx$ 50kHz	7,5 c)

a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT.

b) The carrier shall be modulated using a 50% duty cycle square wave signal.

c) r.m.s., before modulation is applied.

## Important Instructions Before Use

1. **WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
2. **WARNING:** PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of Wrist-type Fully Automatic Digital Blood Pressure Monitor, including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could result.
3. The software identifier refer to the software evaluation report , and the file code is **JYRJ201012001**.

### 4.verify manometer pressure accuracy:

In the power down state, press and hold the "  " button, and theninstall the batteries. Until the LCD screen is full, release the "  " button.

When the LCD screen displays the double zero, the bloodpressure meter is in static state. At this point, 500ml gas capacity, calibrated standard pressure gauge and manual pressure device can be connected to the sphygmomanometer through the sleeve interface of the sphygmomanometer, and manual pressure can be applied to the effective display range of the sphygmomanometer, and then the difference between the reading of the sphygmomanometer and that of the standard pressure gauge can be compared.This mode can be used to verify manometer pressure accuracy.

### 5.Contraindications:

Product is not intended for infants or individuals who cannot express their intentions.

### 6.Intended Use

The digital blood pressure monitor are reusable for clinical and home use and are non-invasive blood pressure measurement systems designed to measure the systolic and diastolic blood pressure and pulse rate of adolescents and adults individual by using a non-invasive technique, which is a well-known technique in the market called the "oscillometric method". it can measure the systolic blood pressure, diastolic blood pressure and pulse rated on up-arm, and the device is reusable for clinical or home use.

7.The patient is the operator:

the PATIENT is an intended OPERATOR.

the PATIENT Do not carry out other maintenance operations except to replace the battery.

8.WARNING:

Do not modify this equipment without authorization of the manufacturer.

9. ESSENTIAL PERFORMANCE Maintenance advice:

Pressure calibration will be carried out when this product leaves the factory. Patients can use the method described in the section "Verify Manometer Pressure Accuracy" to verify the accuracy. If the accuracy deviation is large, please contact the manufacturer to recalibration.

10.Mechanical strength and resistance to heatThe resistance to heat will be retained by device during the EXPECTED SERVICE LIFE of the ME EQUIPMENT.

11.Do not place the blood pressure monitor and cuff at will. It will cause asphyxiation if the child swallows or twine around his neck.

12.The cuff and the case of the blood pressure monitor have been tested for biocompatibility and do not contain allergenic or harmful materials.Please stop using it if allergy occurs during use.

13.Warning:

Non-professionals do not modify the equipment, otherwise it will make the equipment measurement is not accurate.

14.Warning:

Do not expose the equipment for a long time, otherwise it will reduce the performance of the equipment.

15.Warning:

This device is not used for children and pets

16.Clean:

The equipment can be cleaned by lay operator according to the cleaning procedures in the instructions

17.Warning:

Do not use a damaged cuff for blood pressure measurement.

18.Warning:

When measuring with the cuff, if the tester feels seriously uncomfortable, press the button of the blood pressure monitor to deflate the cuff, or remove the cuff directly from the arm.

19. Warning:

If an unexpected reading occurs, the operator can take several more measurements and consult a doctor.

20. Warning:

This equipment is used outside the specified environment, may damage the equipment, and may be inaccurate measurement.

**Correct Disposal of This Product**  
(Waste Electrical & Electronic Equipment)

This marking shown on the product indicates that it should not be disposed with other household waste at the end of its life. To prevent potential harm to the environment or to human health, please separate this product from other types of wastes and recycle it responsibly. When disposing this type of product, contact the retailer where product was purchased or contact your local government office for details regarding how this item can be disposed in an environmentally safe recycling center.

Business users should contact their supplier and check the terms and conditions of the purchasing agreement. This product should not be mixed with other commercial wastes for disposal. This product is free of hazardous materials.

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