

GADNIC

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
NB-1100


Índice

Español	5
Português	31
English	57

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Compresor Nebulizador

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

User's manual

Nebulizer Compressor

Please read this manual thoroughly to learn about the features of this product. Please retain this manual for future reference. Any use of this product other than its intended use is not recommended.

Thank you for purchasing the Compressor Nebulizer.

Intended Use/Indications for Use:

The Compressor Nebulizer system include a compressor and nebulizer. The electrically powered compressor provides compressed air to the supplied pneumatic nebulizer to aerosolize drugs for inhalation by the patient. The nebulizer is driven by the integral air compressor. The device may be used with pediatric (ages 2 years and above) and adult patients in the home, hospital, and sub-acute care settings.

Contraindication:

Please not use the nebulizer device if not follow doctor's advice;
Pentamidine is not applicable to the device.

Principle:

When the compressed air generated by the compression pump is sprayed from the nozzle, the liquid medicine is sucked up through the negative pressure generated between the nozzle and the air tube. The liquid impinges on the upper battle and become a very fine mist and spray out to the outside.

Explanation of the symbol

59

	Manufacturer
	Date of manufacture
	Batch Code
	Unique device identifier
	Storage and Transportation Temperature Limit: -20 °C ~55 °C (-4 °F ~131 °F)
	Type BF Equipment
	Refer to instruction manual/booklet
	Atmospheric pressure limitation
	Storage and Transportation Humidity limitation: 5%~93%RH
	Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products.
IP21	The first num.2: Protected against solid foreign objects of 12,5 mm Ø and greater. The second num.1: Protection against vertically falling water drops
	Caution
	Warning

To ensure the correct use of the product, read all instructions carefully before using the device. Failure to read and observe all precautions could result in personal injury or equipment damage.

CAUTIONS

(OPERATING)

1. The product should never be left unattended when plugged in.
2. The patient is an intended operator, close supervision is necessary when this product is used by, or near children or disabled individuals.
3. Operate the device only as intended. Do not use the device for any other purpose.
4. Not use the unit outdoors or in very dusty environment.
5. Always ensure that the air vents are not blocked or obstructed in any way.
6. Do not disassemble the nebulizer or attempt to repair it; any attempt for repair should be completed only by qualified technical personnel.
7. If your pet's hair or a bug falls into the accessories, please have it cleaned and disinfected.
8. Keep the nebulizer away from your pets.

(ELECTRICAL)

9. Always turn off and unplug the unit after use.
10. When this product is used near other electronic equipment, pay attention to the possible electromagnetic interference between them. The electromagnetic interference may cause the loss or reduction of the basic performance of the machine, and they may need to be separated.

(CLEANING AND STORAGE)

11. Make sure the air filter is clean. If the air filter has discoloration or has been used for more than 70 applications, replace it with a new one.
12. Do not use abrasive cleaners.
13. Do not wrap the power cord around the compressor.
14. Store the device and the components in a clean, safe location.

 WARNING

(OPERATING)

1. Please seek medical attention immediately, if you have any allergic reaction.
2. The nebulizer can work with the solution only, but not with suspensions, emulsions or high viscosity drugs.
3. Operate the device only as intended. Do not use the nebulizer for any other purpose or in a manner inconsistent with these instructions.
4. For type, dose, and regime of medication follow the instructions of your physician or licensed healthcare practitioner.
5. Never use any liquid in the nebulizer other than those prescribed by your doctor. Liquids such as cough medicines or essential oils can harm both the machine and the patient.
6. Do not immerse the compressor in liquid and do not use while bathing. If the unit falls into water, don't touch the device unless it is unplugged, otherwise there is a risk of electric shock.
7. Do not use the unit if it has been dropped, exposed to extreme temperatures or high humidity or damaged in any way.
8. Keep the device and kids accessories out of reach of unsupervised infants and children. The device may contain small accessories that may pose a choking hazard.
9. Do not use in anaesthetic or ventilator breathing circuits.
10. Never use while sleeping or drowsy.
11. Not suitable for use in the presence of flammable anesthetic mixture with air or oxygen or nitrous oxide.
12. Do not operate the device where oxygen is being administered in a closed environment, i.e. oxygen tent or oxygen hood.
13. Do not allow the nebulized medicine to come in contact with your eyes.
14. Always ensure that the air vents are not blocked or obstructed in any way.
15. Please stop using the device immediately if the nebulizer is not working properly like when it makes unusual sounds, or if you feel pain or discomfort while using.
16. Do not expose the unit to direct sunlight, heated or hot surfaces, humid environments, extreme temperatures, strong static electricity or electromagnetic waves.

 WARNING

(OPERATING)

- 17. Not to position the nebulizer to make it difficult to operate the disconnection device.
- 18. Use of accessories or detachable parts other than those specified by the manufacturer may result in unsafe or degraded performance.
- 19. Please do not connect other parts not recommended by the manufacturer to the atomizer to avoid unnecessary wrong connection.
- 20. Please away from children to prevent strangulation due to cables and hoses.

(ELECTRICAL)

- 21. Do not use the compressor (main unit) or power cord while they are wet.
- 22. Do not use when bathing or with wet hands.
- 23. Do not touch the main unit for other than necessary operation such as turning off the power while nebulizing.
- 24. Do not operate the device with a damaged power cord or plug.
- 25. Unplug the power cord from the electrical outlet before cleaning the device.
- 26. If the power cord is damaged or in other circumstances, and need to replace the power cord, contact the professional personnel of the manufacturer. Do not replace the power cord yourself.
- 27. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- 28. PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Compressor Nebulizer, including cab is specified by the manufacture. Otherwise, degradation of the performance of this equipment could result.


 WARNING

(CLEANING AND STORAGE)

- 29. Disconnect the unit from the electrical outlet before cleaning.
- 30. Never immerse the unit in water to clean as it may damage the unit.
- 31. Do not place or attempt to dry the device, components or any of the nebulizer parts in a microwave oven.


NOTE: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Compressor nebulizer unit

Nebulizer kit

Nebulizer Cup

Nebulizer kit

Mouthpiece

Adult Mask

Child Mask

Air Tube

Other included items

Air Filters Package

Owner's Manual

Standard Accessories and Replace Parts
Air Filter (×5 pcs)

NOTE:

The manufacturer has placed a air filter in the nebulizer before leaving the factory.

A air filter expected use 70 applications or turn into yellow before discard. The old air filter discard according to the requirements of the local authorities and regulation.

Whether you are a medical staff or an ordinary user, you can quickly start without professional training just follow the instructions.

Inspection for use:

1. Please inspect whether the air filter add into the compressor nebulizer unit and the air filter keep normal color and in the range of expected lifetime.
2. Inspecting the power cord is in good condition. If it is damaged, please stop to use.
3. Please check whether the compression pump works properly before use. Press the I/O switch, the indicator light will come on and the sound of the compression pump will be heard if the device can be work normally.

Operate instruction:

1. Plug the power cord into the electrical outlet, make sure the I/O (ON/OFF) switch is in the O (OFF) position. (See Figure 1)

NOTE: Compressor can be worked at 220V, 50 Hz.


Figure 1

2. Add the correct amount of prescribed medication to the medication cup, hold the cup steady.(See Figure 2)
3. Place the baftle on the nozzle.(See Figure 3)
4. Place the cap on the medication cup and twist to close.(See Figure 4)
5. Attach the mask or mouthpiece to the medication cup.(See Figure 5)
6. Attach the air tube to the air tube connector of medication cup and nebulizer, gently rotate the air tube plug to ensure it connect securely.(See Figure 6)

Figure 2


Figure 3

Figure 4

Figure 5


Figure 6

7. Insert the mouthpiece into your mouth, or place the mask over the nose and mouth. (See Figure 7)
8. Turn on the nebulizer by pressing the power switch to the I (ON) position and begin treatment.(See Figure 8)


Figure 7


Figure 8

9. Hold the nebulizer kit straightly. Exhale normally through the mouthpiece or mask.
- If treatment needs to be interrupted, simply press power switch to O (OFF) position.
10. When the treatment is complete, turn off the compressor by pressing the power switch to O (OFF) position.
11. Unplug the power cord from the electrical outlet.

⚠ WARNING

For type, dose, and regime of medication, follow the instructions of your physician or licensed healthcare practitioner.

Hold the medication cup correctly. Do not tilt the medication cup so the angle of the medication cup is greater than 45 degrees.

If the nebulizer has any fault during use, unplug it immediately.

The maximum temperature of the atomization chamber is 5 °C higher than the ambient temperature.

⚠ CAUTION

If the unit has been stored at temperatures below 0°C, leave it in a warm place for about 15 minutes before using it.

If the unit has been stored at temperatures above 40°C, leave it in a cool place for about 15 minutes before using it.

Replacing the Air Filter:

To replace the air filter, follow these steps(See Figure 9):

1. Change the air filter only when the unit is off, unplugged and completely cool.
2. Gently pull the air filter cover out of the nebulizer.
3. Remove the air filter from the air filter cover.
4. Place a new air filter in the air filter cover.
5. Return the air filter cover to its place, ensuring that it is firmly in place.

Figure 9

NOTE: Make sure the air filter is clean. If the air filter has changed color, replace it with a new one.

WARNING

Carefully wash your hands before cleaning and disinfecting the nebulizer parts. To reduce the risk of infection, illness, or injury due to contamination, please clean and disinfect all parts of the nebulizer after each treatment. And nebulizer cups are for single patient use only. This will decrease the risk of contamination.

NOTE: Cleaning and disinfection have no effect on the accessories performance. It is advisable to replace the Nebulizer Kit after about 100 treatments. And if them has discoloration or deformation, replace them specified by the manufacturer.

Cleaning

Cleaning the Nebulizer Kit:

1. Disassemble the nebulizer kit into single parts before cleaning;
2. Prepare cleaning solution for the mask, nebulizer cup, mouthpiece, and airway in sequence according to the instructions for using neutral detergent.
3. Submerge the parts in the cleaning solution for 5 minutes. After soaking, use a lint-free cloth to wipe twice. During each wipe, all areas that can be touched need to be wiped and thoroughly cleaned.
4. Completed step C, rinse the parts thoroughly with tap water for 1.5 minutes.
5. Shake off excess water and dry with a lint-free cloth, then place in a clean place to dry for at least 30min.

Enclose of Compressor unit:

1. With power switch in the O (OFF) position, unplug power cord from outlet.
2. Wipe daily with a soft cloth moistened with water or mild, neutral detergent.
3. Do not use any powdered cleaners or soap pads, which may damage the unit and Air Tube.

Disinfecting

Always disinfect the parts before using the product for the first time and after the product has not been used for an extended period of time as well as after the last treatment of the day.

Disinfection by Alcohol::

Follow the instructions provided by the disinfectant manufacturer.

1. Clean according to the instructions before disinfection;
2. Submerge the nebulizer cup and mask in Alcohol solution at room temperature ($23^{\circ}\text{C} \pm 2^{\circ}\text{C}$) for 15 minutes.
3. Remove the parts after disinfection, rinse them thoroughly with tap water. Then dry with a lint-free cloth and place in a clean place to dry for at least 30min.

NOTE: Never clean with benzene, thinner or a flammable chemical.

Disinfect by boiling:

At room temperature, boil the nebulizer cup and mouthpiece parts for 15 to 20 minutes. After boiling, carefully remove the parts, then shake off excess water and allow to air dry in a clean environment.

Please select the disinfection method from the table below:

○ : applicable ×: not applicable

Part	Nebulizer Cup	Mouthpiece	Mask	Air Tube
Disinfection Method	Medium level disinfection	High level disinfection	Medium level disinfection	Cleaning
Boiling	○	○	×	×
Alcohol	○	×	○	×

WARNING

Do not use a microwave oven, dish dryer or hair dryer to dry parts.

Cleaning the compressor nebulizer unit:

1. Making sure the nebulizer compressor is power off.
2. Wipe daily with a soft cloth moistened with water or mild, neutral detergent.
3. Wipe with the soft cloth moistened with water over 5 times if you use neutral detergent in step 2.
4. Air dry and save for next time.

Inspection and maintenance within device shelf life:

1. Check whether should be replace air filter, the replace method please see the chapter “Replacing the Air Filter”.
2. Check the power cord whether in good condition, please stop use if the power cord broken or exposed internal wires.
3. Clean nebulizer compressor unit shell if necessary.
4. Do not open or disassemble the compressor at any time. This action will void the warranty. Repairs service should be carried out only by us authorized technician.

NOTE:

- Do not store the device in extreme hot or cold temperatures, high humidity or under direct sunlight.
- Do not crease or fold the air tube.
- Do not subject the device or the parts to any strong shocks such as dropping on the floor.
- The compressor nebulizer unit is used for same patient or multiple patient.
- Dispose of the device and any used accessories or optional parts according to applicable local regulations.
- Never clean with benzene, thinner or a flammable chemical.
- Do not use the device where it may be exposed to flammable gas.

⚠ Electric shock hazard.

Do not remove compressor cabinet. All disassembly and maintenance must be done by authorized service center..

PROBLEM	CAUSE	SOLUTION
No power to unit when the power switch is on.	The power cord is not plugged into an electrical outlet.	Turn the power switch off. Plug the power cord into an electrical outlet. Turn the device on.
	The fuse was broken	Find professional technicians for repair
No nebulization or low nebulization rate when the power is on.	No medication in the medication cup. Too much or too little medication in the medication cup.	Add the correct amount of prescribed medication to the medication cup.
	The nebulizer kit is not correctly assembled.	Make sure the nebulizer kit is correctly assembled and the inhalation accessory is correctly attached.
	The nozzle is blocked.	Clean the medication cup to remove the blockage.
	The medication cup is tilted at an incorrect angle.	Hold the medication cup correctly. Do not tilt the medication cup so the angle of the medication cup is greater than 45 degrees.
	The air tube is incorrectly attached.	Make sure the air tube is correctly attached to the nebulizer and the medication cup.
	The air tube is folded or damaged. The air tube is blocked.	Make sure the air tube is not folded, kinked or bent. Inspect the air tube for any damage. Replace the air tube if damaged.
The device is abnormally loud.	The air filter cover is incorrectly attached.	Attach the air filter cover correctly. Make sure the air filter cover is not blocked.
The device is very hot.	The compressor is covered.	Do not cover the compressor with any type of cover during use.
	The equipment operates at temperatures above 40°C	Please operate the compressor at normal temperature

SPECIFICSTIONS	
Type:	Compressor Nebulizer
Power supply:	AC 220V, 50 Hz
Input Power:	≤ 80VA
Operating Mode:	Non-Continuous (20min. On, 30min. OFF)
Sound Level:	≤ 65dB(A)
Free Flow Rate:	≥ 6.5L/min
Dimensions:	Approx. 10.5" (w) × 6.2" (h) × 4.1" (d) (267.9 mm × 157.3 mm × 104.6 mm)
Weight:	Approx. 1.55kg (Compressor only)
Operating Temperature / Humidity / Atmospheric Pressure	+5°C to +40°C (+41°F to +104°F)/15% to 90% RH/ 86 kPa to 106 kPa
Storage and Transport Temperature / Humidity / Atmospheric Pressure	-20°C to +55°C (-4°F to +131°F)/5% to 93% RH/ 70 kPa to 106 kPa
Normal Working Pressure:	50kPa~180kPa
Abnormal Working Pressure:	190kPa~320kPa
Functions:	Power cord storage
Classification:	Class II equipment, Type BF applied part
Contents:	1 Compressor nebulizer unit, 5 Air Filters, 1 Owner's manual
Applied Parts:	Mask, Mouthpiece
Accessible Parts:	Shell of Compressor
Fuses:	5JF 1.6AL 250V

NOTE: Although the shell of the compressor is not used as the applied parts, but as the parts that touches the human body, it needs to be designed according to the requirements of the applied parts.

Manufacturer will provide circuit diagrams, component part lists, descriptions, calibration instructions to assist to service personnel in parts repair.

The following disclosures for nebulizer performance are based upon testing that utilizes adult ventilatory patterns and are likely to be different from those stated for paediatric or infant populations

Test parameters of compressor nebulizer unit with nebulizer kit: fully open	
Medication Cup Capacity	8 mL (cc) max.
Appropriate Medication Quantities	4mL (cc)
Residual Volume*	≤ 1mL
D ₅₀ *	3.0μm ± 25%
Equivalent particle size distribution*	The proportion of particles smaller than 5μm is greater than 60%
Nebulization Rate*	≥ 2.5 mL/min

Test parameters of compressor nebulizer unit with nebulizer kit: fully close	
Medication Cup Capacity	8 mL (cc) max.
Appropriate Medication Quantities	4mL (cc)
Residual Volume*	≤ 1mL
D ₅₀ *	3.0μm ± 25%
Equivalent particle size distribution*	The proportion of particles smaller than 5μm is greater than 65%
Nebulization Rate*	≥ 0.1 mL/min

*D₅₀ = Median Particle Diameter

Tested using malvern laser by JOYTECH Healthcare Co., Ltd. Measured with NaCl(0.9%) at ambient temperature of 23 ± 2 °C and humidity of 53%.

Specifications are subject to change without notice.

Parts and accessories not approved for use do not meet these specifications.

Performance may vary with testing conditions, environment and drugs.

Using a solution, suspension, or emulsion different from that recommended by the manufacturer, in particular, a suspension and/or high-viscosity solution, can alter the technical datas, which can then be different from those disclosed by the manufacturer.

Durable periods are as follows, provided the product is used to nebulize 2 ml-5ml of medication 2 times a day for 15 minutes each time at room temperature (23°C). Clean the compressor after each use.

Durable period may vary depending on usage environment.

Compressor Nebulizer expected shelf life: 3 years

Air filter: Turn into yellow before discard

Frequent use of the product may shorten the durable period.

This product complies with the following electromagnetic compatibility standard:
EN 60601-1 Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance

EN 60601-1-2 General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests

EN 60601-1-11 Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

EN ISO 27427 Anaesthetic and respiratory equipment - Nebulizing systems and components.

Correct Disposal of This Product (Waste Electrical & Electronic Equipment)

This marking shown on the product indicates that it should not be disposed with other household waste at the end of its life. To prevent potential harm to the environment or to human health, please separate this product from other types of wastes and recycle it responsibly.

When disposing this type of product, follow your local municipality guidelines for disposal. The Controller, AC Supply and Connection Cord can be disposed or recycled with other electronic equipment. You can contact the retailer where product was purchased or contact your local government office for details regarding how this item can be disposed in an environmentally safe recycling center.

Business users should contact their supplier and check the terms and conditions of the purchasing agreement. This product should not be mixed with other commercial wastes for disposal. This product is free of hazardous materials.

Please report to the manufacturer and the competent authority of the Member State in which you are established about any serious incident that has occurred in relation to this device.

There are no hazardous materials used in the construction of this unit, therefore, there are no known risks associated with disposal of the equipment.

The compressor nebulizer is guaranteed from the date of purchase. If the device does not function properly due to defective components or poor workmanship, we will find a possible solution. This warranty does not cover damage due to improper handling in any way.

The single used nebulizer kit (air tube, nebulizer cup, adult mask, child mask, mouthpiece) is not within the range of warranty, but you can contact retailer or dealer to replace a new nebulizer kit if the package is damaged at the first use.

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1

Guidance and declaration of manufacturer-electromagnetic emissions	
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.	
Emissions test	Compliance
Radiated emission CISPR 11	Group 1
Conducted emission CISPR 11	class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies

Table 2

Guidance and declaration of manufacturer-electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
IMMUNITY test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air
Electrostatic transient/burst IEC 61000-4-4	Power supply line: ± 2 kV input/output lines: ± 1 kV	Power supply line: ± 2 kV
Surge IEC 61000-4-5	line(s) to line(s): ± 1 kV line(s) to earth: ± 2 kV 100kHz repetition frequency	line(s) to line(s): ± 1 kV
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % 25/30 cycles Single phase: at 0° 0 % 300 cycle	0 % 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % 25/30 cycles Single phase: at 0° 0 % 300 cycle
RATED power frequency magnetic fields IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz

Table 2 (continued)

Guidance and declaration of manufacturer-electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
IMMUNITY test	IEC 60601 test level	Compliance level
Conducted disturbances induced by RF fields 61000-4-6	3 V rms 0,15 MHz – 80 MHz 6 V rms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V rms 0,15 MHz – 80 MHz 6 V rms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF EM fields 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

Table 3

Recommended separation distances between portable and mobile RF communications equipment and the device				
The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter			
	80 MHz to 800 MHz $d=1.2\sqrt{P}$	150 MHz to 80 MHz $d=0.6\sqrt{P}$	80MHz to 800MHz $d=0.35\sqrt{P}$	800 MHz to 2.7 GHz $d=0.7\sqrt{P}$
0.01	0.12	0.06	0.035	0.07
0.1	0.38	0.19	0.11	0.22
1	1.2	0.6	0.35	0.7
10	3.8	1.9	1.1	2.21
100	12	6	3.5	7
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.				
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

Table 4

Guidance and declaration of manufacturer-electromagnetic emission							
The device is intended for use in the electromagnetic environment specified below. Thecustomer or the user of the device should assure that it is used in such an environment.							
Proximity fields from RF wireless communications equipment 61000-4-3	Test frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (A/m)
	385	380 to 390	TETRA 400	Pulse modulation b) 18 Hz	1.8	0.3	27
	450	430 to 470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 to 787	LTE Band 13,17	Pulse modulation b) 217 Hz	0.2	0.3	9
	745						
	780						
	810	800 to 960	GSM 800/900,TETRA 800,iDEN 820,CDMA 850,LTE Band 5	Pulse modulation b) 18Hz	2	0.3	28
	870						
	930						
	1720	1700 to 1990	GSM 1800; CDMA 1900; GSM1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28
	5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
	5500						
5785							

If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

Table 4(Continued)

- a) For some services, only the uplink frequencies are included.
 b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
 c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50 % duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distance for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E=6/d\sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

Table 5

Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields			
Proximity magnetic fields IEC 61000-4-39	Test frequency	Modulation	IMMUNITY TEST LEVEL(A/m)
	30 kHz a)	CW	8
	134,2kHz	Pulse Modulation b) 2,1kHz	65c)
	13,56MHz	Pulse Modulation b) 50kHz	7,5 c)
a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT. b) The carrier shall be modulated using a 50% duty cycle square wave signal. c) r.m.s., before modulation is applied.			

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