GADNIC

BALANZA DIGITAL

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BALANZ13

ÍNDICE

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IMPORTANT INFORMATION

IMPORTANT:

Please read this user manual carefully to fully familiarize yourself with your new product before first use.

DEAR CUSTOMER:

Thank you for purchasing this product. To ensure your safety and long-term trouble-free use, please read the following safety instructions carefully:

- Read and carefully keep all instructions.
- Heed all warnings.
- Follow all instructions carefully.
- Do not use this device near water.
- Clean only with a dry cloth.
- Do not install near heat sources, such as radiators, heating regulators, kitchens, and other heat-generating appliances (including amplifiers).

IMPORTANT: COMMERCIAL/PROFESSIONAL USE PROHIBITED

BATTERIES

Operates with 2 AAA batteries of 1.5V each. (Not included)

Battery Replacement

Open the compartment at the back of the scale.

Install the new AAA batteries, please observe the polarities.

OPERATION

SET/Setup Button: Power on, select menu, select unit and confirm.

Down Button

Up Button

INSTRUCTIONS

WEIGHING BODY WEIGHT ONLY

Step 1 - Place the scale on a firm, flat surface and stand on it. Do not use the scale on a carpet as this may affect the accuracy of the measurement.

Step 2 - Proceed as indicated in this manual.

- a) When the LCD screen displays 0.0, the scale will be ready to weigh.

0.0 kg

- b) Step gently onto the glass platform of the scale.
Try to move as little as possible until the LCD screen displays the weight reading, at which point the body weight measurement will be complete.
- c) After reading your weight, step off the scale. The device will automatically power off after a few seconds.

MEASUREMENT: WEIGHT, BODY FAT (%), WATER (%), MUSCLE MASS (%), BONE MASS (KG), BMI, AND CALORIES

To measure your body fat (%), water (%), muscle mass (%), bone mass (kg), BMI, and calories, configure the scale as indicated below:

- When the LCD screen displays "0.0 kg", press the **(SET)** cbutton (setup) for a few seconds to set up your personal data.
Press **(▲)**(up), or **(▼)**(down) to select a number from 0 to 9 as your personal code when the numeric code blinks on the LCD screen, and once you reach your number, press **(SET)** to confirm.
- Press **(▲)** o **(▼)** to select the male or female option when the gender symbol blinks on the LCD screen, and then press the **(SET)** button to confirm.


28

- Press  or  to select your height (80-220 cm) when the number and unit of measurement (cm) blink on the LCD screen, then press the  button to confirm.



- Press  or  to select your age (6-100) when the option (age) and number blink on the LCD screen, then press the  button to confirm. Configuration will finish after this step.

- Step gently onto the scale and keep your bare feet in contact with the electrodes. Stay still to obtain your measurements and you will see your weight, body fat (%), water (%), muscle mass (%), bone mass (kg), BMI, and calories on the LCD screen in a few seconds. The results will repeat three times and then the scale will automatically turn off if not in use.

- Your measurements will be stored under your numeric code if it is set up. Simply confirm your code when using the scale again and the device will display your previous information, so there will be no need to reprogram it.

BODY FAT RANGE IN PERCENTAGE: 5%-85%

Standard body fat ranges for adults

Age	Femenino				Masculino			
	Low fat	Healthy	Overweight	Obesity	Low fat	Healthy	Overweight	Obesity
6-39	<21	21-33	33-39	>39	<8	8-19	19-25	>25
40-59	<23	23-35	35-40	>40	<11	11-22	22-28	>28
≥60	<24	24-36	36-42	>42	<13	13-25	25-30	>30

BODY HYDRATION RANGE IN PERCENTAGE: 5%-85%

Comparison table between measured hydration and standard hydration level

Age	Female			Male		
	Low	Normal	Hihg	Low	Nonnat	Hihg
6-15	<57	57-67	>67	<58	58-72	>72
16-30	<47	47-57	>57	<53	53-67	>67
31-60	<42	42-52	>52	<47	47-61	>62
≥60	<37	37-47	>47	<42	42-56	>56

MUSCLE MASS (%): 13% - 80%

Female: <34% - Male: >40%

BONE MASS (KG): 0.5-8.0 KG

	Female			Male		
Weight	<45 kg	45-60 kg	>60 kg	<60 kg	60-75 kg	>75 kg
	>1.8 kg	>2.2 kg	>2.5 kg	>2.5 kg	>2.9 kg	>3.2 kg

IMPORTANT INFORMATION

ERROR INDICATIONS

Overload indication:

The letters "Err" will appear on the LCD screen when the weight placed on the scale exceeds 180 kg. Remove the weight from the scale to avoid damage.

Low battery indication:

The letters "LO" (low) will appear on the LCD screen. Open the battery compartment cover and replace the existing batteries with new ones, correctly placing the negative and positive polarities.

Failed measurement:

Failed measurement is indicated by "Err!".

SPECIFICATIONS

- High precision strain gauge system.
- Capacity: 5 kg - 180 kg (unit: kg).
- Attention: This scale will turn on when 6 kg or more are placed on it. After it is turned on, the minimum weight this scale can weigh is 5 kg. Precision: 100 g.
- Overload indication: (Err).
- Low battery indication: (Lo).
- LCD backlight.
- Memory for 8 users (storage of personal information).
- Measurements: weight, body fat (%), water (%), muscle mass (%), bone mass (kg), BMI, and calories.

ATENCIÓN

- Please remove shoes and socks and keep bare feet in contact with the electrodes of the scale to measure your weight, body fat (%), water (%), muscle mass (%), bone mass (kg), BMI, and calories.
- If the scale does not function correctly, check the batteries. Replace them with new batteries if necessary.
- Use soft paper towels with alcohol or glass cleaner to clean the surface of the scale if it is dirty. Do not use soap or other chemicals. Keep the scale away from water and extreme temperatures, whether high or very low.
- The scale is a high precision measuring device. Do not jump or stomp on the scale with excessive force, and avoid disassembling it. Handle with care to avoid damage when moving.
- This scale is for household use only and is not suitable for professional use. Measurements of body fat (%), water (%), muscle mass (%), bone mass (kg), BMI, and calories are for reference only. Consult a doctor before starting any diet or exercise program.
- Caution: Scale may become slippery when wet! The glass platform can become quite slippery when wet, so be sure the glass and your feet are dry before use. Never stand on the scale with wet feet.

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