MASAJEADOR OCULAR LED

LED EYE MASSAGER · MASSAJADOR
OCULAR LED

LEDMAS56

ANTES DE USAR, APLIQUE AGUA O GEL SOBRE EL PARCHO · APPLY WATER OR GEL TO THE PATCH
BEFORE USE · ANTES DE USAR, APLIQUE ÁGUA OU GEL NO PATCH

LEA ATENTAMENTE LAS INSTRUCCIONES ANTES DE SU USO · READ CAREFULLY BEFORE USE · LEIA
ATENTAMENTE ANTES DE UTILIZAR

01

ÍNDICE • CONTENTS • ÍNDICE

ESPAÑOL	03
ENGLISH	05
PORTUGUÊS	07

ESPAÑOL

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MASAJEADOR OCULAR LED

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USER MANUAL

LED EYE BEAUTY MASSAGER

BEFORE USE, APPLY WATER OR GEL TO THE WORKING PATCH OR EYE BAGS TO ENSURE THE BEST EFFECT.

READ THE INSTRUCTIONS CAREFULLY BEFORE USE. KEEP THIS MANUAL FOR FUTURE REFERENCE.

PRODUCT SPECIFICATIONS

PRODUCT NAME	Eye beauty massager
MATERIAL	ABS
BATTERY CAPACITY	90 mAh
CHARGING TIME	30 minutes
USAGE TIME	90 minutes
INCLUDED ACCESSORIES	Eye massager, charging cable, user manual

OPERATING MODES

BUTTON PRESS	FUNCTION
Long press	Power on / off
Short press 1st time	ES fixed frequency + red light
Short press 2nd time	ES variable frequency + red light
Short press 3rd time	ES variable frequency only
Short press 4th time	Red light only

PREPARATION AND USE

1. INITIAL CHECK

After long-pressing the button, check that it enters working mode. If not, charge for **1.5-2 hours** before use.

2. CLEAN THE FACE

Remove makeup thoroughly, especially around the eye bag area.

3. APPLY WATER OR GEL

After cleaning, apply water or gel to the working patch or eye bags to ensure the ES current can conduct to the skin. This step is essential for optimal results.

4. USAGE NOTES

- The product is not waterproof; do not rinse with water. If cleaning is needed, use a damp cloth or towel.
- Default mode on startup is ES fixed frequency + red light, with ES intensity set to level 1.
- Every time the mode changes, intensity resets to level 1.
- Adjust ES intensity with short press (3 levels: low, medium, high).
- Automatic timer: automatically shuts off after **10 minutes**.

5. CHARGING INDICATORS

- Red light: charging in progress.
- Green light: fully charged. Recommend using a **5V / 1A** charger.

- Low battery reminder: red light flashes **3 times**.

FREQUENTLY ASKED QUESTIONS

- **Can't feel the massage or only one side is massaging?** Apply water or gel to the patch or skin. Press the instrument firmly against the skin to ensure the patch adheres tightly to the skin.
 - **How often should I use it?** It can be used every day, with each session recommended for **10 minutes**.
-

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