GADNIC

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PASSIVE LEG EXERCISER

INSTRUCTION MANUAL

READ THE INSTRUCTIONS CAREFULLY BEFORE USE. KEEP THE MANUAL FOR FUTURE REFERENCE.

WELCOME

We invite you to sit and walk with this new fitness device. The main dynamic is very simple: while your legs move, your heart moves too, helping to keep it healthy. Thanks to it, you can combat the effects of reduced mobility and sedentary lifestyles comfortably and effectively, without leaving home. After regular use, you will begin to notice benefits throughout the body, especially in the legs and joints, promoting blood circulation, muscle activation and oxygenation. Its use is highly recommended and has no contraindications.


PRODUCT FEATURES

PRODUCT TYPE	Tabletop passive leg exerciser, for the feet
SPEED LEVELS	30
AUTOMATIC MODES	3 (P1, P2, P3)
MANUAL MODE	1
REMOTE CONTROL	Included
MAGNETS	On the foot support surfaces
FRONT LEGS	Height adjustable
AUDIO INPUT	USB, with built-in speakers
WEIGHT	4.2 kg
DIMENSIONS	20 x 41 x 32.7 cm

CONTROL PANEL

SPEED DISPLAY	TIME DISPLAY
ON/OFF function	Switching on and off
"P" function	Default modes
Hand function	Manual mode
Speed control	Decreases speed (works in Manual mode)
Speed control	Increases speed (works in Manual mode)

PRODUCT PARTS

- 1 Pedals with magnetic surface and support magnets
- 2 Control panel with LED display
- 3 Foot securing straps
- 4 Height-adjustable front legs
- 5 USB audio input / power cable

PRECAUTIONS FOR USE

- Do not unplug the device with wet hands, or while it is switched on.
- Do not place a heavy load on the power cable and do not crush it; do not move the device by pulling the cable, and do not wrap the cable around the device.
- Do not use a damaged power cable or switch.
- Do not spill any liquid on the equipment or use it in damp areas; do not touch it with wet hands during operation.
- Do not let children play with the equipment. It must not be used by children or people with disabilities without the necessary supervision of an adult.
- Do not place the machine on a slippery surface.
- In case of fire, do not wet the device if it is switched on.
- Use the equipment in a well-ventilated room.
- This product is suitable for family use; do not use it for professional purposes.
- Check that it is not damaged before use, and stop using it if you notice anything unusual.
- Switch off the device when not in use or once you have finished using it.
- If you feel uncomfortable while using the device, stop immediately and consult your doctor.
- Do not cover the device during operation (e.g. with a blanket or pillow), as it could overheat and cause damage.
- Do not use the device while sleeping.
- This device is designed for use on the feet; do not use it on other parts of the body.
- Do not use the device right before or after eating; wait at least 1 hour.
- Do not use the device if the ambient temperature exceeds 40 °C. If moving it from a cold environment, wait at least one hour before switching it on, as internal condensation may form and affect its operation.
- Do not insert your hands or feet between the pedals and the machine at any time.
- Do not stand on the machine; it is not a toy.
- In the event of a power outage, disconnect the plug from the electrical outlet.
- To unplug the device, pull the plug rather than the cable.
- To prevent fire, damage or electric shock, the correct voltage must be used.
- Do not expose the device to high temperatures; store it in a dry place.
- If you have any health issues, consult your doctor before using the device.
- If you notice the device is not working properly, do not attempt to repair, disassemble or modify it. Unauthorised repairs can cause damage and void the warranty; contact customer service.

POTENTIAL BENEFITS

If used regularly, the following results can be obtained: it helps to reactivate blood flow throughout the body; it helps to prevent and reduce chronic venous ailments; and it helps to improve and prevent cardiovascular health by activating the circulatory system.

GENERAL INFORMATION

Recommended for people with a sedentary lifestyle. Sitting or lying down for long hours can cause blood to circulate unevenly from the legs to the heart. With 3 hours of daily use, the benefits are similar to those achieved by walking: it helps reduce swelling in the legs and ankles, increases the amount of oxygen reaching the muscles, and helps reduce ankle pain.

DAILY USE

It is very useful, comfortable and easy to use. With regular use it helps improve posture in the workplace and reduces postural fatigue syndrome, helps eliminate lower back pain and leg heaviness, and its continued use can help prevent venous discomfort.

REMOTE CONTROL

BUTTON	FUNCTION
1	Ignition
2	Starting up
3 / 7	Speed up
4	Mode change (P1-P2-P3, once)
5 / 6	Speed down
8	Mode change
9	Speed level 10
10	Speed level 30
11	Speed level 20
12	Pause movement

OPERATING MODES

MODE	DESCRIPTION
Manual	"HR" will appear on the LED indicator
Automatic	P1, P2 and P3, selectable
P1	Speed varies every minute between levels 1 and 10, looping up and down
P2	Speed varies every minute between levels 11 and 20, looping up and down
P3	Speed varies every minute between levels 21 and 30, looping up and down

THE IMPORTANCE OF WALKING

Walking regularly, more than 3 hours a week, can help reduce the risk of cardiovascular disease: it strengthens the heart, arteries, muscles and bones of the legs and feet; supports digestion and metabolism; helps clear the mind and reduce fatigue; and helps lower blood pressure, cholesterol and abdominal fat.

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