GADNIC

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MASAJ028

ÍNDICE

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“User Manual”.

Product Features

- 6 Massage Heads.
- 6 Massage Speeds.
- 6 massage modes.
- Extend the electrical line (3 meters).
- High temperature output.
- Massage heads can massage all parts of the body.
- Practical to eliminate fatigue and relax.
- Elongated massager for longer reach.

Mode of use

- Insert the plug into the socket.
- Press “ ” to turn on the unit.
- Press the “+”, “-” button to adjust the massage speed.
- Press the “M” button to adjust the massage mode.
- Press “ ” to activate the heating mode.
- Every program will be deactivated after 15 minutes of use.
- Pull out the plug properly after use.

Caution

- Discontinue use when the machine is very hot.
- Disconnect the equipment in a cautious way, without pulling the cable.
- Pay attention to children or handicapped persons who use it, they must be supervised for the use of the equipment.

Massage head instruction

	Round massage head: the whole body part can be massaged.
	Heated massage head: the heat from the massager will help improve fatigue and tension in the working area.
	6-node massage head: deep massage for back and waist.
	Massage head with essential oil: apply essential oil or cream on your body and massage with the massager.
	Foot massage head: pleasant head for the feet and additional pain relief.
	Acupuncture head: stimulates particular acupuncture points.

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